

The Highest Log Reduction Available in a Hand Sanitizer!

Yet, goes on light and leaves employees hands feeling soft and silky.

Delivered as an atomized spray for maximum coverage

71% Ethanol formula

Dual moisturizers help keep hands soft and healthy even with repeated use

Meets all FDA (GRAS) and USDA requirements for food handling

NSF listed E3 classification

Kosher, Pareve and Halal certified

Wide variety of dispensing options











Atomized spray saturates fingernails, cuticles, creases and crevices of the fingers and hands where pathogens hide.





See product label for handling and directions for use. See SDS for additional information.

NSF Registration

Nonfood Compounds Program Listed (E3) Registration #137057.

Active Ingredients

Ethyl Alcohol 71.0000 % Inactive ingredients: Water, isopropanol, glycerin, glycerol monolaurate

Warning: Flammable **Color:** Clear Liquid

Pounds Per Gallon: 7.25 (0.84kg/L)

Other information

Store in a cool, dry place below 104° F (40° C)

Uses

- Sanitizes hands to help reduce bacteria that potentially can cause disease
- Helps to prevent cross-contamination
- Helps prevent drying of skin
- Recommended for repeated use

Directions

- Wet hands thoroughly with spray, covering entire surface of hands, especially the fingernails and cuticles
- Rub hands together until dry



Areas Most Often Missed with Inadequate Hand Washing.

The graphic on the left shows the areas most often missed with inadequate hand washing. Alpet E3 Plus Hand Sanitizer Spray is the most effective way to kill germs where they hide: under fingernails, under cuticles, and in the creases and crevices of fingers and hands.

Efficacy

Acinetobacter baumannii ATCC 19606
Bacteroides fragilis ATCC 25285
Candida albicans ATCC 10231
Candida glabrata ATCC 26512
Enterobacter aerogenes ATCC 13048
Enterococcus faecalis ATCC 29212
Enterococcus faecium ATCC 19434
Escherichia coli ATCC 25922
Escherichia coli 0157:H7 ATCC 43895

Haemophilus influenzae ATCC 19418 Klebsiella pneumoniae ATCC 4352 Klebsiella oxytoca ATCC 43165 Listeria monocytogenes ATCC 7644 Micrococcus luteus ATCC 7468 Proteus mirabilis ATCC 7002 Pseudomonas aeruginosa ATCC 15442 Pseudomonas aeruginosa ATCC 27853 Serratia marcescens ATCC 14756

Staphylococcus aureus ATCC 29213 Staphylococcus aureus MRSA ATCC 33592 Staphylococcus epidermidis ATCC 12228 Staphylococcus haemolyticus ATCC 29970 Staphylococcus hominis ATCC 27844 Staphylococcus saprophyticus ATCC 15305 Streptococcus pneumoniae ATCC 6302 Streptococcus pyogenes ATCC 19615

*Alpet E3 Plus was tested against Salmonella typhi using a chlorine equivalency test and shows efficacy superior to a 200 ppm chlorine solution.

Configurations

SA10036 48x4 oz. finger-pump bottles SA10012 6x1000 mL countertop pump bottles SA10039 6x1000 mL VersaClenz cartridges

SA10014 4x1 gallon bottles

SA10013 50 gallon drum (bulk for refill)

SA20002 500mL Secondary Container

SA20001 1000 mL Secondary Container

SA20000 One Gallon Secondary Container

USP20092 Secondary Container Tag (Dark Blue)







*Salmonella typhi ATCC 6539



Dispensers

AD10048 VersaClenz Touchless Dispenser (White - single)

AD10048C VersaClenz Touchless Dispenser (White - case of 6) AD10048B VersaClenz Touchless Dispenser (Black - single)

AD10048BC VersaClenz Touchless Dispenser (Black - case of 6)

MD10030 VersaClenz Manual Dispenser (White - single)

MD10030C VersaClenz Manual Dispenser (White - case of 6)

MD10030B VersaClenz Manual Dispenser (Black - single)

MD10030BC VersaClenz Manual Dispenser (Black - case of 6)

MD10200 Gallon Wall Bracket - Stainless Steel MD10006 1-2 Knockout Wrist-Activated Dispenser - Stainless Steel

MD10105 EZ-Step Wall-Mounted Foot-Activated Dispenser - Stainless Steel

MD10012 EZ-Step Portable Foot-Activated Dispenser - Stainless Steel

LT10008 Laminated Bilingual Instructional Hand Sanitizing Poster















